

MSc Research Practicum

CONFIGURATION MANUAL

**Enhancing Leukemia Diagnosis with Synthetic
Data and Explainable Deep Learning
Architectures**

x23373385

Ibrahim Malik

MSc Data Analytics

Group B

This configuration manual provides all necessary parameters and settings to replicate the experimental setup for the leukemia diagnosis research project.

Table of Contents

1. Environment Requirements.....	3
Platform Setup.....	3
Required Libraries.....	3
Installation Commands.....	3
2. Dataset Configuration.....	4
Dataset Source.....	4
Directory Structure.....	4
Image Preprocessing.....	4
Data Augmentation (Training).....	4
3. Model Configurations.....	5
3.1 GAN Architectures.....	5
DCGAN Configuration.....	5
WGAN Configuration.....	5
cGAN Configuration.....	5
3.2 CNN Architecture.....	5
3.3 Vision Transformer (ViT).....	6
3.4 Hybrid CNN-ViT.....	6
4. XAI Configuration.....	7
4.1 Grad-CAM Settings.....	7
4.2 SHAP Configuration.....	7
4.3 LIME Settings.....	7
5. File Organization.....	8
Code Structure.....	8
6. Training Configuration.....	9
Callbacks (All Models).....	9
Evaluation Metrics.....	9
7. Common Issues & Solutions.....	9
Memory Issues.....	9
Training Issues.....	9
GAN Issues.....	9
Colab Issues.....	9

1. Environment Requirements

Platform Setup

- **Platform:** Google Colab Pro (required for enhanced GPU/TPU access)
- **Python Version:** 3.12
- **Storage:** Google Drive integration for dataset access
- **Hardware:** GPU runtime (Tesla T4)

Required Libraries

All libraries use latest available versions:

- TensorFlow/Keras
- NumPy
- OpenCV-Python
- Matplotlib
- Scikit-learn
- SHAP
- LIME
- Scipy

Installation Commands

```
pip install tensorflow keras numpy opencv-python matplotlib scikit-learn shap lime scipy
```

2. Dataset Configuration

Dataset Source

- **Primary Source:** The Cancer Imaging Archive (TCIA)
- **URL:** [C_NMC_2019 Dataset: ALL Challenge dataset of ISBI 2019 \(C-NMC 2019\) - The Cancer Imaging Archive \(TCIA\) Public Access - Cancer Imaging Archive Wiki](#)
- **Dataset Name:** C-NMC (Children's Hospital of Philadelphia)

Directory Structure

/content/drive/MyDrive/leukemia_detection_project/




Image Preprocessing

- **Target Size:** 224×224 pixels
- **Normalization:** [0,1] scale (divide by 255.0)
- **Color Mode:** RGB (3 channels)
- **Data Type:** float32

Data Augmentation (Training)

- Rotation: ±20 degrees
- Width/Height Shift: ±10%
- Shear: ±10%
- Zoom: ±10%
- Horizontal Flip: Enabled
- Rescale: 1/255

3. Model Configurations

3.1 GAN Architectures

DCGAN Configuration

- **Training:** 225 epochs, batch_size=32
- **Image Normalization:** [-1,1] range
- **Latent Dimension:** 100
- **Generator:** 7 Conv2DTranspose layers (7×7→224×224, tanh output)
- **Discriminator:** 2 convolutional layers (stride=2, sigmoid output)
- **Optimizers:**
 - Discriminator: Adam(lr=0.0002, $\beta_1=0.5$, clipvalue=1.5)
 - Generator: Adam(lr=0.0001, $\beta_1=0.5$)
- **Loss:** Binary cross-entropy

WGAN Configuration

- **Training:** 200 epochs, batch_size=32
- **Latent Dimension:** 500
- **Training Protocol:** 5 critic iterations per generator update
- **Gradient Penalty:** $\lambda_{gp}=10$
- **Generator:** Dense→7×7, then 5 Conv2DTranspose layers (tanh output)
- **Critic:** Identical to DCGAN discriminator (no sigmoid)
- **Optimizer:** Adam(lr=0.0001, $\beta_1=0.5$) for both networks

cGAN Configuration

- **Training:** 250 epochs, batch_size=32
- **Latent Dimension:** 500
- **Classes:** 0=healthy, 1=leukemic
- **Generator:** Dense→7×7, then 5 Conv2DTranspose + concatenated class labels
- **Discriminator:** Receives both image and label inputs
- **Optimizers:**
 - Generator: Adam(lr=0.0002, $\beta_1=0.5$)
 - Discriminator: Adam(lr=0.00005, $\beta_1=0.5$, clipvalue=1.0)

3.2 CNN Architecture

- **Input Shape:** (224, 224, 3)
- **Architecture:**
 1. Conv2D(32, 3×3) + ReLU → MaxPooling2D(2×2)
 2. Conv2D(64, 3×3) + ReLU → MaxPooling2D(2×2)
 3. Conv2D(128, 3×3) + ReLU → MaxPooling2D(2×2)
 4. Flatten → BatchNormalization
 5. Dense(128) + ReLU → Dropout(0.5)
 6. Dense(1) + Sigmoid
- **Optimizer:** Adam(lr=0.0001)
- **Loss:** Binary cross-entropy
- **Training:** 50 epochs, batch_size=32

3.3 Vision Transformer (ViT)

- **Input Shape:** (224, 224, 3)
- **Patch Size:** 16×16 (P=16)
- **Number of Patches:** 196 ($N=(224/16)^2$)
- **Embedding Dimension:** 256 (D=256)
- **Transformer Layers:** 4 (L=4)
- **Attention Heads:** 4 (h=4)
- **Head Dimension:** 64 (dk=64)
- **MLP Hidden Sizes:** [512, 256]
- **Activation:** GELU
- **Optimizer:** Adam(lr=1e-4)
- **Loss:** Categorical cross-entropy
- **Training:** 30 epochs, batch_size=32

3.4 Hybrid CNN-ViT

- **CNN Branch:** 4 conv blocks (filters: 32, 64, 128, 256)
- **ViT Branch:** Processes CNN feature maps, 16×16 patches -> 9 patches
- **Patch Dimensions:** 65,536 (16×16×256)
- **Transformer Blocks:** 6 (L=6)
- **Fusion Method:** Element-wise addition
- **Classification Head:** [512, 256] hidden units, GELU activation
- **Training:** 50 epochs, batch_size=32, Adam(lr=1e-4)

4. XAI Configuration

4.1 Grad-CAM Settings

- **Target Layers:**
 - CNN: 'conv2d_2' (last convolutional layer)
 - Hybrid: 'conv2d_3'
- **Class Indices:** 0=Healthy, 1=Leukemia
- **Visualization:** Alpha=0.6, Colormap='jet'
- **Quality Threshold:** heatmap_std > 0.1

4.2 SHAP Configuration

- **Explainer Type:** Partition Explainer (primary), Deep Explainer (fallback)
- **Parameters:**
 - num_samples=50 (background images)
 - max_evals=500 (model evaluations)
 - batch_size=32
- **Masker:** Inpainting (preferred), blur(1,1) (fallback)


4.3 LIME Settings

- **Parameters:**
 - num_samples=1000 (perturbed samples)
 - num_features=5 (top superpixels to highlight)
 - positive_only=True
 - hide_rest=False
- **Segmentation:** quickshift(kernel_size=4, max_dist=200, ratio=0.2)

5. File Organization

Code Structure

Code/



6. Training Configuration

Callbacks (All Models)

- **Early Stopping:** monitor='val_loss', patience=5
- **Model Checkpoint:** save_best_only=True
- **ReduceLROnPlateau:** factor=0.2, patience=3, min_lr=1e-7

Evaluation Metrics

- Accuracy, Precision, Recall, F1-Score, AUC-ROC
- Sensitivity, Specificity
- FID, SSIM, Inception Score (for GANs)

7. Common Issues & Solutions

Memory Issues

- **Symptoms:** ResourceExhaustedError, CUDA out of memory
- **Solutions:** Reduce batch_size to 16 or 8, enable GPU memory growth

Training Issues

- **Symptoms:** Model not converging, NaN values
- **Solutions:** Check learning rates (try 1e-4 or 1e-5), verify data normalization

GAN Issues

- **Symptoms:** Mode collapse, training instability
- **Solutions:** Balance G/D training rates, use label smoothing, try WGAN for stability

Colab Issues

- **Symptoms:** Session timeouts, runtime disconnections
- **Solutions:** Use Google Colab Pro, save checkpoints frequently, remount Drive if needed